

Technical features

- Aluminium frame.
- Adjustable and padded upholstery.
- TOP SPIN 1 and 3 height and tilt adjustable pulling down back-rest.
- TOP SPIN 2 height adjustable back-rest could be pulled down (only if higher than 30 cm).
- Tilt, height, forward & backward adjustable seat to optimise the set-up (TOP SPIN 1 and 3).
- 24" or 26 quick-release rear wheels equipped with tubulars or pneumatic tyres, and aluminium, rubber coated or titanium pushrims.
- 80 mm front castors.
- 80 mm height adjustable rear castor.
- Height and depth adjustable foot-rest.
- Rear wheel camber: 18°.
- Weight:

Top Spin 1	10,6 kg
Top Spin 2	9,3 kg
Top Spin 3	11,3 kg


 (standard set-up)



Take up less space



Seat with adjustable set-up (TOP SPIN 1 and 3)



TOP SPIN 1 AND 3 STANDARD SET-UP TRACE

Order form

- TOP SPIN 1 TOP SPIN 2 TOP SPIN 3
- Frame length (TOP SPIN 1 only): **(A)** 58 cm 64 cm (TOP SPIN 2 and 3 only 51 cm)
- Seat width: 36 cm 38 cm 40 cm 42 cm
- Seat depth: 38 cm 40 cm 42 cm
- Front seat/floor height: **(B)** from 34 cm to 54 cm: cm
- Rear seat/floor height: **(C)** from 30 cm to 40 cm: cm
- Back-rest height:cm
- Back-rest: pulling down fixed
- Rear tyres: tubulars pneumatic tyres tubeless tyres
- Distance between seat and foot-rest: **(D)**cm
- Upholstery colour: black cornflower blue red
- Frame colour: (see "Frame colours" leaflets)
- Rear wheels with aluminium pushrims: 24" 26"
- Distance between hub centre and back-rest (from 8 cm to 18 cm): **(E)**cm

Options

- Spinergy Sport wheels: 24" 25" 26"
- Polycarbonate spoke guards: with **OFFCARR** logo without **OFFCARR** logo
- Foam polyurethane lined cushion (75 Shore), h: 3cm 5cm
- PVC leg keeper belt: length cm and widthcm
- Rubber coated pushrims: Titanium pushrims:
- Front grip handle: detachable fixed
- right left



Via dell' Artigianato II, 29
 35010 Villa del Conte (PD) Italy
 Tel. +39 049 9325733 - Fax +39 049 9325734
 www.offcarr.com - offcarr@offcarr.com



Dealer address

Date: _____
 Order number: _____
 Customer ref.: _____
 Notes: _____
