

# Tennis


# Top Spin

## Technical features


- Aluminium frame.
- Adjustable and padded nylon back-rest upholstery.
- TOP SPIN 1, 3 and 4 height adjustable back-rest.
- TOP SPIN 1 and 3 pulling down back-rest (only if higher than 30 cm).
- Tilt, height, forward & backward adjustable seat to optimise the set-up (TOP SPIN 1 and 3).
- 24", 25" and 26" quick-release rear wheels, with aluminium, rubber coated or titanium pushrims.
- 80 mm front castors.
- 80 mm height adjustable rear castor.
- Height and depth adjustable foot-rest.
- Rear wheel camber: 18°.
- Weight: Top Spin 1 10,6 kg  
Top Spin 4 10 kg  
Top Spin 3 11,3 kg.  
(standard set up)



Adjustable seat  
(TOP SPIN 1 and 3)



TOP SPIN 4 STANDARD SET-UP TRACE



TOP SPIN 1 AND 3 STANDARD SET-UP TRACE

## Order form

- TOP SPIN 1     TOP SPIN 3     TOP SPIN 4  
 -Frame length (TOP SPIN 1 only): (A)  58 cm     64 cm (TOP SPIN 4 and 3 only 51 cm)  
 -Seat width:  36 cm     38 cm     40 cm     42 cm  
 -Seat depth:  38 cm     40 cm     42 cm  
 -Front seat/floor height: (B) from 34 cm to 54 cm: ..... cm  
 -Rear seat/floor height: (C) from 30 cm to 40 cm: ..... cm  
 -Back-rest height: .....cm  
 -Back-rest:  pulling down     fixed     tilt adjustable  
 (Top spin 1 and 3)  
 -Distance between seat and foot-rest: (D) .....cm  
 -Upholstery colour:  black     cornflower blue     red  
 -Frame colour: .....  
 (see "Frame colours" leaflets)  
 Rear wheels with aluminium pushrims:  24"     25"     26"  
 Distance between hub centre and back-rest  
 from 8 cm to 18 cm: (E) .....cm (compulsory indication for TOP SPIN 4)  
 Anti-tip device:  single     twin

## Options

- Spinergy Sport wheels:  24"     25"     26"  
 Polycarbonate spoke guards:  with **OFFCARR** logo     without **OFFCARR** logo  
 Foam polyurethane lined cushion (75 Shore), h:  3cm     5cm  
 Leg keeper belt: lenght ..... cm and width ..... cm  
 Rubber coated pushrims :  Titanium pushrims:   
 Front grip handle:  detachable     fixed  
 right     left



Via dell' Artigianato II, 29  
35010 Villa del Conte (PD) Italy  
Tel. +39 049 9325733 - Fax +39 049 9325734  
www.offcarr.com - offcarr@offcarr.com

Dealer address

Date: \_\_\_\_\_  
 Order number: \_\_\_\_\_  
 Customer ref.: \_\_\_\_\_  
 Notes: \_\_\_\_\_