

Dealer \_\_\_\_\_

### Key


standard

option

- L** Seat width 34 36 38 40 42 44 46  
**P** Seat depth 36 38 40 42 44  
**HS** Back-rest height 27 29 31 33 35 37 39  
**HA** Front seat/floor height 46 48 50  
**HP** Rear seat/floor height 39 41 43 45 47  
**DP** Foot-rest/seat distance 31 33 35 37 39 41 43 45 47 49

cm

## VENUS ELITE



### Frame colour

- Shiny titanium  Pearl white  Glossy black  Matt black  
 Natural titanium  
 ..... (see "Colours" catalogue)

①

### Upholstery

②

- Tension adjustable nylon back-rest  
 Tension adj. breathable back-rest  
 Without back-rest

### Single tubular titanium foot plate

③

### Pushing handles

④

- Open  Carbon fibre closed  Tilt adjustable Carbon fibre closed

- Without any pushing handles  Fixed pushing handles  Height adjustable pushing handles  Foldable pushing handles

### Distance between the hub centre and back-rest

⑤

### Back-rest tilt

⑥

### Frame length

⑦

### Seat with shock absorbers

⑧

- 4 cm  9 cm  
 5 cm  11 cm  
 7 cm (standard)

- 89°  
 91° (standard)  
 93°

- 42,5 (std)  
 46,5  
 .....

cm

User's weight : ..... kg  
 It is not possible to have two anti-tip devices

### Side guards

⑨

### "L" detachable arm-rest - h 22 cm

⑩

### Brakes

⑪

### Revolving anti-tip device

⑫

- Aluminium with carbon fibre wing  Detachable  
 Aluminium with straight carbon fibre  Detachable

- Tip-up  
**Elbow rest**  
 Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)





- Push brake  Pull brake  Scissor brake

- Right  Left \_\_\_\_\_

### Notes

Dealer

<p><b>Rear wheels</b> (12)</p> <ul style="list-style-type: none"><li><input type="checkbox"/> 24" Light (24 spokes)</li><li><input type="checkbox"/> 24" Star 6x3 (18 spokes)</li><li><input type="checkbox"/> 24" (36 spokes)</li><li><input type="checkbox"/> 24" Mountain-bike w/Citybike tyres</li><li><input type="checkbox"/> 24" Mountain-bike w/Offroad tyres</li><li><input type="radio"/> 24" Ultralight (24 spokes) black</li><li><input type="radio"/> 24" Spinergy</li><li><input type="radio"/> 25" Spinergy</li><li><input type="radio"/> 25" Light (24 spokes) black</li></ul> <p><b>Tyres</b></p> <ul style="list-style-type: none"><li><input type="checkbox"/> Pneum. Light Run 24"</li><li><input type="checkbox"/> Pneum. Marathon Plus</li><li><input type="checkbox"/> Pneum. high press. smooth</li><li><input type="checkbox"/> Pneum. high press. plugged</li><li><input type="radio"/> Solid 24" (not available with Ultralight)</li></ul>	<p><b>Handrims</b> (13)</p> <ul style="list-style-type: none"><li><input type="checkbox"/> Aluminium silver</li><li><input type="checkbox"/> Aluminium dark</li><li><input type="checkbox"/> Chrome plated steel 24"</li><li><input type="radio"/> Titanium</li><li><input type="radio"/> Rubber coated</li><li><input type="radio"/> Tween Grepp 24"</li><li><input type="radio"/> Max Grepp 24"</li><li><input type="radio"/> Max Grepp Ergo 24" (**)</li><li><input type="radio"/> Silicon cover 24"</li></ul> <p>(**) not suitable for wheels Star 6x3 and Mountain bike</p>	<p><b>Spoke guards</b> (14)</p> <ul style="list-style-type: none"><li><input type="radio"/> Transparent</li><li><input type="radio"/> White</li><li><input type="radio"/> Black</li><li><input type="radio"/> Multicolour</li><li><input type="radio"/> Black plumed red</li><li><input type="radio"/> Black plumed blue</li><li><input type="radio"/> Black plumed green</li><li><input type="radio"/> Alphabet</li></ul> <p><b>Rear wheel camber</b> (15)</p> <ul style="list-style-type: none"><li><input type="checkbox"/> 1°</li><li><input type="checkbox"/> 2°</li><li><input type="checkbox"/> 3°</li><li><input type="checkbox"/> 4°</li></ul>	<p><b>Front castors</b> (16)</p> <ul style="list-style-type: none"><li><input type="checkbox"/> 80 (3") solid</li><li><input type="radio"/> 80 (3") solid with lights</li><li><input type="checkbox"/> 100 (4"x1.1/4) solid with plastic rim</li><li><input type="radio"/> 100 (4"x1) solid with lights</li><li><input type="checkbox"/> 100 (4"x1) solid with aluminium rim</li><li><input type="checkbox"/> 100 (4"x1.1/4) solid with aluminium rim</li><li><input type="checkbox"/> 125 (5"x1.1/4) solid with plastic rim</li><li><input type="checkbox"/> 125 (5"x1) solid with aluminium rim</li><li><input type="checkbox"/> 125 (5"x1.1/4) solid with aluminium rim</li></ul>
---	--	---	--

<p><b>Pockets</b> (17)</p>  <ul style="list-style-type: none"><li><input type="checkbox"/> Small under the seat</li></ul>  <ul style="list-style-type: none"><li><input type="radio"/> Big under the seat</li></ul>  <ul style="list-style-type: none"><li><input type="checkbox"/> On the back-rest</li></ul>	<p><b>Lined cushion in foam polyurethane</b> (18)</p>  <ul style="list-style-type: none"><li><input type="radio"/> 3 cm high</li><li><input type="radio"/> 5 cm high _____</li></ul> <p><b>Anodized set</b> (19) (Rear wheel rims, rear wheel hubs, forks)</p> <ul style="list-style-type: none"><li><input type="checkbox"/> Silver</li><li><input type="checkbox"/> Black</li><li><input type="checkbox"/> Lime Green</li><li><input type="checkbox"/> Light blue</li><li><input type="checkbox"/> Red</li><li><input type="checkbox"/> Orange</li><li><input type="checkbox"/> Cyclamen light purple</li></ul> <p>(only 24" Light wheels)</p>	<p><b>Options</b></p> <ul style="list-style-type: none"><li>20 <input type="checkbox"/> "M" velcro strip (hook) on seat</li><li>21 <input type="checkbox"/> "F" velcro strip (loop) on seat</li><li>22 <input type="radio"/> Padded rigid seat h 5</li><li>23 <input type="radio"/> Chest harness</li><li>24 <input type="radio"/> Chest belt</li><li>25 <input type="radio"/> Abductor belt</li><li>26 <input type="radio"/> Catheter holder pocket</li><li>27 <input type="radio"/> Rx <input type="radio"/> Lx Crutch holder</li><li>28 <input type="radio"/> Pump for pneumatics</li><li>29 <input type="radio"/> Suitcase holder</li><li>30 <input type="radio"/> 45° pelvic belt</li><li>31 <input type="radio"/> 45° padded pelvic belt</li><li>32 <input type="radio"/> Black holdall</li><li>34 <input type="radio"/> Umbrella holder</li><li>35 <input type="radio"/> Front protective covers</li></ul>
--	---	---

**Notes**

---

---

---

---

---

---

---