


Dealer \_\_\_\_\_

**Key**       standard       option

# FENICE

- L** Seat width       34    36    38    40    42    44
- P** Seat depth       36    38    40    42    44
- HS** Back-rest height       27    29    31    33    35    37    39
- HA** Front seat/floor height       46    48    50
- HP** Rear seat/floor height       40    42    44    46
- DP** Foot-rest/seat distance       35    37    39    41    43    45    47    49      cm



- Colore telaio**      ①
- Metallic grey       Pearl white       Metallic light blue
  - Glossy black       Cyclamen light purple
  - .....(see "Colours" catalogue)

- Upholstery**      ②
- Tension adjustable nylon back-rest
  - Tension adj. breathable back-rest
  - Without back-rest

- Front frame**      **Single tubular foot-plate**      ③
- V
  - Y
  - Closed in carbon
  - Open (height adjustable +/- 2 cm only with Y shaped front frame)

- Pushing handles**      ④
- Without any pushing handle
  - Standard pushing handle
  - Height adj. pushing handles
  - Foldable pushing handles

- Distance between the hub centre and back-rest**      ⑤
- 
- 7 cm       13 cm
  - 9 cm (standard)       15 cm
  - 11 cm

- Back-rest tilt**      ⑥
- 
- 89°
  - 91° (standard)
  - 93°

- Front angle**      ⑦
- 
- 90° (standard) (A = 16 cm)
  - 95° (A = 22 cm)
- (measures are approximate for a standard arrangement)

- Rear wheel camber**      ⑧
- 1°       2°
  - 3°       4°

- Side guards**      ⑨
- 
- Carbon fibre       Plastic
  - € 80,00

- "L" arm-rest, detachable, h 22cm      ⑩
- 
- Elbow-rest**
- Tubular
  - Straight (26 cm long)
  - Classic (23 cm long)
  - Classic (30 cm long)

- Brakes**      ⑪
- 
- Standard
  - Scissor brakes

- Revolving anti tip device**      ⑫
- 
- Right
  - Left

Notes \_\_\_\_\_

