

Dealer _____


Key

standard

option

ZODIAC

- L** Seat width 34 36 38 40 42 44
- P** Seat depth 36 38 40 42 44
- HS** Back-rest height 27 29 31 33 35 37 39
- HA** Front seat/floor height 46 48 50
- HP** Rear seat/floor height 40 42 44 46
- DP** Foot-rest/seat distance 35 37 39 41 43 45 47 49 cm



Frame colour

- Metallic grey Pearl white Metallic light blue
- Glossy black Cyclamen light purple
- (see "Colours" catalogue)

①


Upholstery

②

- Tension adjustable nylon back-rest
- Tension adj. breathable back-rest
- Without back-rest

Single tubular foot-plate

③




(not height adjustable)

- Open Closed

Pushing handles


④



- Without pushing handles Standard pushing handles Height adjustable pushing handles

Distance between the hub centre and back-rest


⑤



- 7 cm 13 cm
- 9 cm (standard) 15 cm
- 11 cm

Back-rest tilt


⑥



- 89°
- 91° (standard)
- 93°

Front angle

⑦



- 83° (standard) (A=16cm)
- 75° (A=23cm)


Rear wheel camber

⑧

- 1° 2°
- 3° 4°

Side guards


⑨



- Aluminium straight Aluminium with carbon fibre wing - detachable
- Aluminium with carbon fibre wing

"L" arm-rest, detachable, h 22cm


⑩



- Tip-up
- Elbow-rest**
- Tubular
- Straight (26 cm long)
- Classic (23 cm long)
- Classic (30 cm long)

Brakes


⑪



- Standard
- Scissor

Revolving anti-tip device

⑫



- Right
- Left

Notes

Dealer _____

Rear wheels ⑫

- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes)
- 24" (36 spokes)
- 24" Mountain-bike with Citybike tyre
- 24" Mountain-bike with Offroad tyre
- 24" Ultralight (24 spokes) black
- 24" Spinergey
- 25" Spinergey
- 25" Light (24 spokes) black

Tyres

- Pneum. Light Runv 24"
- Pneum. Marathon Plus
- Pneum. high press. smooth
- Pneum. high press. plugged
- Solid 24" (not available with Ultralight)

Handrims ⑬

- Aluminium silver
- Aluminium dark
- Chrome plated steel 24"
- Titanium
- Rubber coated
- Tween Grepp 24"
- Max Grepp 24"
- Max Grepp Ergo 24" (**)
- Silicon covers 24"

(**) not available with Star 6x3 and Mountainbike wheels

Spoke guards ⑭


- Transparent
- White
- Black
- Multicolour
- Black plumed red
- Black plumed blue
- Black plumed green
- Alphabet

Front castors ⑮


- 80 (3"x1) solid
- 80 (3"x1) solid with lights
- 100 (4"x1.1/4) solid with plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1) solid with aluminium rim
- 100 (4"x1.1/4) solid with aluminium rim
- 125 (5"x1.1/4) solid with plastic rim
- 125 (5"x1) solid with aluminium rim
- 125 (4"x1.1/4) solid with aluminium rim

Pockets ⑯


Small under the seat




Big under the seat



On the back-rest



Lined cushion in foam polyurethane ⑰



- 3 cm high
- 5 cm high

Anodized set ⑱
(Rear wheel rims, rear wheel hubs, forks)

- Black
- Silver
- Lime Green
- Light blue
- Red
- Orange
- Cyclamen light purple

(only 24" Light wheels)

Options

- 19 "M" velcro strip (hook) on seat
- 20 "F" velcro strip (loop) on seat
- 21 Padded rigid cushion h 5
- 22 Chest harness
- 23 Chest belt
- 24 Abductor belt
- 25 Catheter holder pocket
- 26 Rx Lx Crutch holder
- 27 Rx Lx Double crutch holder
- 27 Pump for pneumatics
- 28 Suitcase holder
- 29 45° pelvic belt
- 30 45° padded pelvic belt
- 31 Black holdall
- 33 Umbrella holder
- 34 Front protective covers

Notes
