


Dealer \_\_\_\_\_

**Key**     standard     option

# ALHENA

<b>L</b> Seat width	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46
<b>P</b> Seat depth	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46
<b>HS</b> Back-rest height	<input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41
<b>HA</b> Front seat/floor height	<input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50
<b>HP</b> Rear seat/floor height	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46
<b>DP</b> Foot-rest/seat distance	<input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47

cm



<b>Frame colour</b> <span style="float: right;">①</span> <input type="checkbox"/> Shiny titanium <input type="checkbox"/> Glossy black <input type="checkbox"/> Pearl white <input type="checkbox"/> ..... <small>(see "Colours" catalogue)</small>	<b>Upholstery</b> <span style="float: right;">②</span> <input type="checkbox"/> Tension adjustable nylon back-rest <input type="checkbox"/> Tens.adjustable breathable backrest <input type="radio"/> Without backrest <input type="checkbox"/> Without seat	<b>Upholstery colour</b> <span style="float: right;">③</span> <input type="checkbox"/> Black <input type="radio"/> Coloured..... <small>(available colours: orange, red, green, cornflower blue, light blue)</small> <input type="radio"/> Jacquard..... <small>(specify pattern)</small> <input type="radio"/> Jacquard..... + Ipad Sleeve+ Pochette <small>(specify pattern)</small>
---	--	---

**Not detachable foot-rest** ④  
(with DP=43 cm and P=40 cm)  
**(standard)**

**Not detachable foot-rest**  
(with DP=43 cm and P=40 cm)

**Detachable foot-rest**

(approx.measurements)

**Front Frame** ⑤

STANDARD

Y  
(only in case of not detachable foot-rest)

**Closed carbon single foot-plate** ⑥

Foot-holding belts

**Carbon Automatic closure foot-plate in carbon** ⑦

Foot-holding belts

**Separate aluminium profile tilt adjustable foot-plates** ⑧

forward     backward

Foot-holding belts

**Single aluminium profile tilt adjustable foot-plate** ⑨

forward     backward

turnable & liftable upwards aside

Foot-holding belts

**Side guards** ⑩

Aluminium with carbon fibre wing

Plastic

Curved aluminium

Straight aluminium

**Tip-up, detachable sport arm-rest h. 22 cm** ⑪

**Elbow-rest**

Tubular  
 Classic (23 cm long)

**Tip-up height adjustable from 22 to 28 cm "L" arm-rest** ⑫  
(detachable, open with sideguards)

Not detachable

**Elbow-rest**

Tubular  
 Straight (26 cm long)  
 Classic (23 cm long)  
 Classic (30 cm long)

**Rear wheels** ⑬

24" Ultralight (24 spokes) black  
 24" Light (24 spokes)  
 24" Star 6x3 (18 spokes)  
 24" (36 spokes)  
 24" Mountain-bike CITYBIKE tyre  
 24" Mountain-bike OFFROAD tyre  
 24" Spinerger  
 25" Spinerger  
 25" Light (24 spokes) black

**Tyres**

Pneumatic Marathon Plus  
 Pneum. high press. plugged 24"  
 Pneum. high press. smooth  
 Solid 24"(not avail.with Ultralight)

**Handrims** ⑭

Aluminium silver  
 Aluminium dark  
 Chromium-plated steel 24"  
 Rubber coated  
 Tween Grepp 24"  
 Max Grepp 24"  
 Max Grepp Ergo 24"<sup>(\*)</sup>  
 Titanium  
 Silicon covers 24"

(\*) not suitable for wheels Star 6x3 and Mountain Bike

Dealer

**Spoke guards** (15)

- Transparent
- White
- Black
- Multicolour
- Black plumed red
- Black plumed blue
- Black plumed green
- Alphabet

**Front castors** (16)

- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) solid plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1) solid with alum. rim
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) solid plastic rim
- 125 (5"x1) solid with alum. rim
- 125 (5"x1.1/4) solid with alum. rim
- 150 (6") pneumatic
- 150 (6") solid

**Revolving anti-tip device** (17)

- ☆ Right
- ☆ Left




**Unbalancing Pedal**

- ☆ Right
- ☆ Left



**Extended wheel plate** (18)




(Not available with scissor brakes)


**Back-rest pushing handles** (19)

- Standard  Without

- Height adjustable



**Distance between hub centre and back-rest** (20)



- 3,5 cm
  - 6 cm (standard)**
  - 8,5 cm
- (with HP=44 cm)

**Anodized set** (21)

(Rear wheel rims, rear wheel hubs, forks and fittings)

- Black
- Silver
- Lime Green
- Light blue
- Red
- Orange
- Cyclamen light purple


(24" Light wheels only)

**Lined cushion in foam polyurethane** (22)




- 3 cm high
- 5 cm high

**Pockets** (24)



- Small under the seat



- On the back-rest

**Wheel camber** (25)

- 1°
- 2°
- 3°
- 4°

**Brakes** (23)

- Classic - short handle
- Classic - pull to brake short handle
- Classic - long handle
- Scissor

**Notes**

**Options**

- 26  "M" velcro strip (hook) on seat
- 27  "F" velcro strip (loop) on seat
- 28  Seat/back-rest connection flap
- 29  Padded rigid cushion h 5
- 30  Chest harness
- 31  Chest belt
- 32  Abductor belt
- 33  Catheter holder pocket
- 35  Back-rest stretched bar
- 36  Back-rest anatomic board
- 37  Rx  Lx Crutch hold./single
- 37  Rx  Lx Crutch hold./double
- 38  Pump for pneumatics
- 39  Suitcase holder
- 40  45° pelvic belt
- 41  45° padded pelvic belt
- 42  Black holdall
- 44  Umbrella holder

**9,7 kg configuration**

- Wheelchair 38x40x31 cm
- Nylon upholstery
- Frame colour: shiny titanium
- Automatic closure foot-rest
- Fixed back-rest pushing handles
- Aluminium sideguards with carbon fibre wings
- Scissor brakes
- Ultralight wheels
- Pneumatic high pressure smooth tyres
- Ultralight handrims
- 100 mm front castors
- Without pockets
- Max load: 100 kg