

Dealer _____



Key
 standard

 option


HAPPY STANDING

User measurements ①

- A - Armpit height cm from the ground
 B - Waist height cm from the ground
 C - Beneath kneecap height cm from the ground
 D - User height cm from the ground
 E - Palm height cm from the ground


Thoracic support height ②

- h 25 cm
 h 20 cm



	24"	26"	28"	32"
User height D	80-100	90-110	110-130	125-140
Frame inclination	5°/30°	5°/30°	5°/30°	5°/30°
Rear wheel camber	4°	4°	4°	4°
Max user load (kg)	50 kg	50 kg	50 kg	75 kg
Total width L	58	64	64	70
Total length P	80	98	98	98

approx measurements, unit: cm

Big wheels ③

- 24"
 26"
 28"
 32" _____

Frame colour ④

- Pearl red
 Orange
 Pastel yellow
 Pearl white
 Metallic grey

Upholstery colour ⑤

- Blue
 Yellow
 Orange
 Lime green
 Red

Handrims ⑥

- Aluminium
(24" silver or dark;
26" - 28" - 32" dark)
 Rubber coated
(24" - 26" - 28")
 Chromium-plated steel
(24")
 Max Grepp
(24")





Spoke guards ⑦

- Transparent
 Black
 White
 Multicolour
 Black plumed blue
 Black plumed green
 Black plumed red
 Alphabet _____

Available only with 24" wheels

Options ⑧

- Polycarbonate table
 with fixed fastening
 Polycarbonate table
 with pantograph connection

 Abductor

 Abductor with pelvis support

 Anatomic
 abductor

 Seat

 Seat with pelvis support
