

Dealer \_\_\_\_\_

Key  standard  option

**ARYA** (seat depth adjustable only)  **ARYA PLUS** (seat width and depth adjustable)

**L** Seat width 20 22 24 26  
24 26 28 30 32  
28 30 32 34 36

**P** Seat depth 20 22 24 26 28  
30 32 34 36 38  
 (adjustable +2, +4, +6 cm from any chosen measurement)


**HS** Back-rest height 25 27 29 31 33 35 37 39 41

**HA** Front seat/floor height 40 42 44 46 48

**HP** Rear-seat/floor distance 36 38 40 42 44 46

**DP** Foot-rest/seat distance ..... cm (min. 15 cm, max 35 cm for HA=40 cm)  
 (min. 15 cm, max 41 cm for HA=46 cm) cm

# ARYA



### Frame colour

- Pearl red  Pearl white  Metallic pink  
 Metallic light blue  Orange FLUO  Green FLUO  
 ..... (see "Colours" catalogue)  
 Two-colour version a) PEARL WHITE sitting frame  
 b) colour of the wheel adjustment bar: ..... (see "Colours" catalogue)

①

### Upholstery

- Tension adjustable nylon  
 Tension adj. breath.back  
 Without back-rest

②

### Upholstery colour


(Nylon only)

- Black  Light blue  
 Red  Green  
 Conflower blue  Orange

③

### Front frame

④




### Foot-rest

⑤

- Not detachable foot-rest  
 Detachable foot-rest

### Single aluminium foot-rest

⑥




- Foot-holding belts

### Separate aluminium foot-rest

Only for L>32 cm (standard front frame)  
 Only for L>28 cm (abducted front frame)



- Foot-holding belts




Open

Liftable single height & depth adjustable foot-plate

Single choice of footplate for front open frame

- Standard pushing handles  Straight height adjustable pushing handles  
 Abducted height adjustable pushing handles  Height adjustable pushing bar



⑧

### Side-guards

- Plastic  
 Carbon fibre  
 Aluminium



⑨

### Height adjustable arm-rest

Elbow-rest


- Classic (23 cm long)  
 Classic (30 cm long)  
 Straight (23 cm long)



⑩

### Distance between the hub centre and back-rest


⑪



- 3 cm (extended wheel plate)  
 0 cm  
 3 cm (standard)  
 5 cm

### Back-rest tilt


⑫



- 89°  
 91° (standard)  
 93°

### Revolving anti-tip device

⑬



- Right  
 Left  
 Right & Left

### Note

Dealer

Mod. OFARYA+R1-EN Send to : OFF CARR S.R.L. Via Dell'Artigianato II, 29 - 35010 Villa del Conte (Padova) Italy Tel: +39 049 9325733 Fax: +39 049 9325734 e-mail: offcarr@offcarr.com - www.offcarr.com

**Rear wheels** (14)

24" Light (24 spokes) (\*)  
 22" Light (24 spokes)  
 20" Light (20 spokes)

**Tyres**

Pneum. high pressure smooth  
 Pneum. high pressure plugged  
 Pneum. plugged (x1"3/8)  
 Marathon Plus 22" or 24"  
 Solid

(\*) Only with P=36 cm and HP from 40 to 46 cm

**Rear wheel camber** (18)

1°    3°    5°

(\*) Available only with 24" or 22" wheels

**Handrims** (15)

Aluminium silver  
 Chrome plated steel  
 Rubber Coated  
 Tween Greep (\*)  
 Max Grepp  
 Max Grepp Ergo  
 Silicon covers

**Spoke guards** (16)


Trasparent if 22"  
 White  
 Black  
 Multicolour  
 Black plumed red  
 Black plumed blue  
 Black plumed green  
 Alphabet

**Front castors** (17)

80 (3") solid with lights  
 80 (3") solid  
 100 (4") solid with lights  
 100 (4") solid  
 100 (4") solid with alum. rim  
 125 (5") solid  
 125 (5") solid with alum.rim  
 150 (6") pneumatic  
 150 (6") solid


**Head-rest** (19)

**Padded**




Fixed  
 Height adjustable

**Moulded**



Fixed  
 Height adjustable


**Orbital**



Moulded    Anatomic

Pls specify "h" & "d" where requested  
Standard: h = 20 cm   d = 15 cm  
h = ..... cm   d = ..... cm


**Abductor** (20)



1    3    5    7 cm    ..... cm  
 Pls choose or specify height from seat

**Lateral supports** (21)

Swing Away




Right  
 Left

Height adjustable (Pls specify h)

**Table** (22)

Trasparent

Pls specify the fastening  
 single    double




(only with arm-rest)

**Anodized set** (23)  
(Rear wheel rims&hubs and forks)

Black  
 Silver  
 Lime Green  
 Light blue  
 Red  
 Orange  
 Cyclamen light purple

**Lined cushion** (24)  
in foam polyurethane



3 cm high  
 5 cm high

**Brakes** (25)

Classic (inside the sideguard)  
 Pull to brake short handle inside the sideguard  
 Drum  
 Drum + classic (inside the sideguard)

**Options**

26  "M" velcro strip (hook) on seat  
 27  "F" velcro strip (loop) on seat  
 28  Padded rigid seat  
 29  Chest harness  
 30  Chest belt  
 31  Abductor belt  
 32  Pump for pneumatics  
 33  Back anatomic board  
 34  45° pelvic belt  
 35  45° padded pelvic belt  
 36  Black holdall  
 38  Umbrella holder  
 39  Front protetive covers

**Note**

-----

-----

-----

-----

-----

-----

-----

-----

-----

-----

Prices are not binding and may change without notice  
Prices refer to 2016 price-list