


Dealer

Key standard option

ALHENA

L Seat width	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
P Seat depth	<input type="checkbox"/> 34 <input type="checkbox"/> 36 <input type="checkbox"/> 38 <input checked="" type="checkbox"/> 40 <input type="checkbox"/> 42 <input type="checkbox"/> 44 <input type="checkbox"/> 46	
HS Back-rest height	<input type="checkbox"/> 29 <input type="checkbox"/> 31 <input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41	
HA Front seat/floor height	<input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input checked="" type="checkbox"/> 50	
HP Rear seat/floor height	<input type="checkbox"/> 38 <input type="checkbox"/> 40 <input type="checkbox"/> 42 <input checked="" type="checkbox"/> 44 <input type="checkbox"/> 46	
DP Foot-rest/seat distance	<input type="checkbox"/> 33 <input type="checkbox"/> 35 <input type="checkbox"/> 37 <input type="checkbox"/> 39 <input type="checkbox"/> 41 <input type="checkbox"/> 43 <input type="checkbox"/> 45 <input type="checkbox"/> 47	cm



Frame colour
 Shiny titanium Glossy black
 Pearl white Natural titanium
 Neon green Neon orange

 (see "Colours" catalogue)

Upholstery
 Tension adjustable nylon back-rest
 Tens.adjustable breathable backrest
 Without backrest
 AMEKH custom postural backrest
(fill up the proper order form)

Upholstery colour
 Black
 Coloured.....
 Available colours:
 NYLON: orange, red, green, cornflower blue, light blue;
 BREATHABLE: dry green, red, grey;

Not detachable foot-rest ④
 (with DP=43 cm and P=40 cm)
 (standard)

Not detachable foot-rest
 (with DP=43 cm and P=40 cm)

Detachable foot-rest
 84° (std)
 76°
(approx.measurements)

Front Frame ⑤
 STANDARD
 Y
(only in case of not detachable foot-rest)
 V
(only in case of not detachable foot-rest & without any automatic closure footplate)

Separate aluminium profile tilt adjustable foot-plates ⑧
 forward backward
 Foot-holding belts

Closed carbon single foot-plate ⑥
 Foot-holding belts

Single aluminium profile tilt adjustable foot-plate ⑨
 forward backward
 up & outward
 Foot-holding belts

Automatic closure foot-plate in carbon ⑦
(note: up to seat width 46cm)
 Foot-holding belts

Unica tubolare aperta in titanio ⑩
 posizionata in avanti posizionata indietro(std)

Side guards ⑩
 Aluminium with carbon fibre wing
 Plastic
 Curved aluminium (same colour or black)
 Straight aluminium (same colour or black)

Tip-up, detachable sport arm-rest h. 22 cm ⑪

Elbow-rest
 Tubular
 Classic (23 cm long)

Tip-up height adjustable from 22 to 28 cm "L" arm-rest ⑫
(detachable, open with sideguards)
 Not detachable

Elbow-rest
 Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)

Rear wheels ⑬
 24" Ultralight (24 spokes) black
 24" Light (24 spokes)
 24" Star 6x3 (18 sps) silver or black
 24" (36 spokes)
 24" Mountain-bike CITYBIKE tyre
 24" Mountain-bike OFFROAD tyre
 24" Spinergy
 25" Spinergy
 25" Light (24 spokes) black
Tyres
 Pneumatic Marathon Plus
 Pneum. high press.plugged 24"
 Pneum. high press. smooth
 Solid 24"(not avail.with Ultralight)

Pushrims ⑭
 Aluminium silver
 Aluminium dark
 Chromium-plated steel 24"
 Rubber coated
 Titanium
 Carbolife CURVE L24**
 CL Curve L TetraGrip 24**
 CL Quadro rubber ctd 24**
 Carbolife GEKKO 24**
 Silicon covers
* not suitable for wheels Star 6x3 and Mountain Bike

