



Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 1/2

2020, January

Dealer

Key standard option

CHILDREN 3000

24* 26* 28* 30 32 34 36

L Seat width

* Not available with detachable foot-rests

CHILDREN 3000 PLUS

L Seat width (adjustable within the selected range)

RANGE 1 24* 26* 28* 30 32 34
 RANGE 2 26* 28* 30 32 34 36 38

* Not available with detachable foot-rests

P Seat depth

28 30 (adjustable from 28 to 34 cm)
 32 34 36 38 (adjustable from 32 to 38 cm)

HS Back-rest height

25 27 29 31 33 35 37 39 41

HA Front seat/floor height

40 42 44 46 48

HP Rear seat/floor height

36 38 40 42 44 46


DP Foot-rest/seat distance

17* 19* 21* 23 25 27 29 31 33
 35 37 39

* Not available with detachable foot-rests

cm

CHILDREN 3000



Frame colour

- Metallic light blue Orange Pearl red Metallic pink
- Neon orange Neon green
- (see "Colour" catalogue)

①

Upholstery

- Tension adjustable nylon back-rest
- Tension adj. breathable b.
- Without back-rest
- Without seat

②

Upholstery colour


(only nylon)

- Black Light blue
- Red Green
- Cornflower blue Orange


③

Front frame

- Standard (Not suggested for L < 32cm)



- Abducted




④

Foot-rest

- Detachable
- Not detachable
- Elevating (detachable) Subject to technical approval

⑤


Separate aluminium profile foot-plates



- Foot-holding belts

⑥

Single aluminium profile foot-plate




- Foot-holding belts

⑦

Side guards


- Plastic (standard)** (Not available with reclining back-rest)



- Aluminium (curved, same colour)

⑧

Height adjustable "T" arm-rest




Elbow-rest


- Classic (23 cm long)
- Classic (30 cm long)
- Straight (26 cm long)

⑨


Standard pushing handles



Straight height adjustable pushing handles



Abducted height adjustable pushing handles



⑩

Rear wheels

- 24" Light (24 spokes) (*)
- 22" Light (24 spokes)
- 20" Light (20 spokes)

Tyres

- Pneum. high pressure smooth
- Pneum. high pressure plugged
- Solid
- Marathon Plus 22" or 24"

(*) Only with P = 34 or 36 cm and HP from 40 to 46 cm

⑪

Handrims

- Aluminium silver
- Chromium-plated steel
- Dark aluminium
- Rubber coated
- Carbolife CURVE L 24"
- Carbolife GEKKO 24"
- Silicon covers

⑫

Spoke guards

- Transparent
- Black
- Alphabet
- Fantasy
- code pattern:
- (see catalogue of spokeguards)

⑬

Front castors

- 80 (3") solid with lights
- 80 (3") solid
- 100 (4") solid with lights
- 100 (4") solid
- 100 (4") solid with alum. rim
- 125 (5") solid
- 125 (5") solid with lights
- 125 (5") solid with alum. rim
- 150 (6") pneumatic
- 150 (6") solid

⑭

Rear wheel camber

- 1° 3° 5°

⑮

Dealer _____

<p>Revolving anti-tip device (16)</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> Right and left</p>	<p>Extended wheel plate (17)</p> <div style="text-align: center;"> </div> <p style="text-align: center;">Compulsory with reclining back-rest</p>	<p>Lined cushion (18) in foam polyurethane</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> 3 cm high <input type="checkbox"/> 5 cm high</p>	<p>Single drive (19) (Only with 0° rear wheel camber) (Only with 22" or 24" wheels)</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> Right <input type="checkbox"/> Left</p>
<p>Head-rest (20) Padded</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> Fixed <input type="checkbox"/> Height adjustable</p> <p>Moulded</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> Fixed <input type="checkbox"/> Height adjustable</p> <p>Orbital</p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <input type="checkbox"/> Moulded </div> <div style="text-align: center;"> <input type="checkbox"/> Anatomic </div> </div> <p style="font-size: 8pt;">Specify "h" & "d" where requested Standard: h = 15 cm d = 10 cm h = cm d = cm</p>	<p>Abductor (21)</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> 1 <input type="checkbox"/> 3 <input type="checkbox"/> 5 <input type="checkbox"/> 7 cm <input type="checkbox"/> cm Choose or specify the height from the seat cover</p>	<p>Lateral supports (22)</p> <div style="text-align: center;"> </div> <p><input type="checkbox"/> Right <input type="checkbox"/> Left Specify the height "h"</p>	<p>Table (23)</p> <p><input type="checkbox"/> Transparent <input type="checkbox"/> Black soft padded</p> <p style="font-size: 8pt;">specify the fastening: <input type="checkbox"/> single <input type="checkbox"/> double</p> <div style="text-align: center;"> (Only with arm-rest) </div>
<p>Anodized set (24) (Rear wheel rim&hub and front fork)</p> <p><input type="checkbox"/> Black <input type="checkbox"/> Silver <input type="checkbox"/> Light blue <input type="checkbox"/> Red <input type="checkbox"/> Orange <input type="checkbox"/> Cyclamen light purple</p>		<p>Brakes (25)</p> <p><input type="checkbox"/> Classic - short handle <input type="checkbox"/> Classic - pull to brake <input type="checkbox"/> Classic - long handle</p> <p><input type="checkbox"/> Right single brake <input type="checkbox"/> Left single brake <input type="checkbox"/> Drum <input type="checkbox"/> Drum + classic</p>	
<p>Reclining back-rest (26)</p> <p>from 90° to 150° (compulsory: - straight or abducted height adjustable - pushing handles - stretched bar - extended wheel plate)</p> <div style="text-align: center;"> </div>		<p>Notes</p> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/> <hr style="border-top: 1px dashed black;"/>	

All prices refer to the price list at time of publishing, are not binding and may be amended without notice.