

Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____ User's weight Kg _____

PAGE 1/2

2020, January


Dealer

Key

standard

option

- L** Seat width 34 36 38 40 42 44 46
P Seat depth 36 38 40 42 44
HS Back-rest height 27 29 31 33 35 37 39
HA Front seat/floor height 46 48 50
HP Rear seat/floor height 39 41 43 45 47
DP Foot-rest/seat distance 31 33 35 37 39 41 43 45 47 49
 cm



Frame colour

- Shiny titanium Glossy black
 Pearl white Matte black
 Natural titanium
 (see "Colours" catalogue)
 Neon green Neon orange

①

Upholstery

- Tension adjustable nylon back-rest
 Tension adj. breathable back-rest
 Without back-rest

②

Upholstery colour

- Black
 Coloured (nylon).....
 colours: orange, red, green, cornflower, light blue
 Coloured (breathable).....
 colours: red, dry green, grey

Single tubular foot-plate

③

- Open Carbon fibre closed
 Aluminium closed
 Tilt adjustable Carbon fibre closed
 Footrest colour: black (standard)
 Foot-rest colour.....


Pushing handles

④

- Without pushing handles Standard pushing handles Height adjustable pushing handles

Distance between the hub centre and back-rest


⑤



- 4 cm
 5 cm
 7 cm (standard)
 9 cm
 11 cm

Back-rest tilt


⑥



- 89° 91° (standard) 93°

Seat with shock absorbers

⑦



User weight : kg
 It is not possible to have two anti-tip devices

Side guards

⑧

- Aluminium with fixed carbon fibre wing Straight aluminium
 Aluminium with detachable carbon fibre wing

Detachable "L" arm-rest - h 22 cm

⑨

- Tip-up Elbow-rest
 Tubular
 Straight (26 cm long)
 Classic (23 cm long)
 Classic (30 cm long)


Brakes

⑩

- Push brake (stand.)
 Pull brake
 Composite brake
 Scissor brake

Revolving anti-tip device

⑪



- Right
 Left

Notes

Order N° _____ Quotation
 Date _____ Signature _____
 Customer Ref. _____

PAGE 2/2
 2020, January

VENUS

Dealer

Rear wheels ⑫

- 24" Ultralight (24 spokes)
- 24" Light (24 spokes)
- 24" Star 6x3 (18 spokes, silver or black)
- 24" (36 spokes)
- 24" Mountain-bike with Citybike tyre
- 24" Mountain-bike with Offroad tyre
- 24" Spinergy
spokes colour.....(black stand.)
- 25" Spinergy
spokes colour.....(black stand.)
- 25" Light (24 spokes) black

Tyres

- Pneum. Marathon Plus
- Pneum. high pressure Smooth
- Pneum. high pressure Plugged 24"
- Pn.Schwalbe Right Run 24"
- Solid 24" not available for Ultralight wheels

Pushrims ⑬

- Aluminium silver
- Aluminium dark
- Chromium-plated steel 24"
- Rubber coated
- Titanium
- CarboLife CURVE L 24" *
- CL Curve L TetraGrip 24" *
- CL Quadro rubber ctd 24" *
- CarboLife GEKKO 24" *
- Silicon covers 24"

* not suitable for wheels Star 6x3 and Mountain bike

Spokes guards ⑭

- Transparent 24"
 - Black 24"
 - Alphabet 24"
 - Fantasy
code pattern
- (see catalogue of spokeguards)


Rear wheel camber ⑮

- 1°
- 2°
- 3°
- 4°

Front castors ⑰


- 80 (3") solid
- 80 (3") solid with lights
- 100 (4"x1.1/4) solid with plastic rim
- 100 (4"x1) solid with lights
- 100 (4"x1.1/4) solid with alum. rim
- 125 (5"x1.1/4) solid with plastic rim
- 125 (5"x1.1/4) solid with alum. rim
- 125 (5"x1) solid with lights

Backrest bar ⑱




- Standard:
P(bar) 7,5cm;
H(bar) 19cm;
dia 19mm.
- Custom:
P.....cm;
H.....cm;
dia.....mm.


Pockets ⑲



Small under the seat



Big under the seat



On the back-rest

Lined cushion
in foam polyurethane ⑳



- 3 cm high
- 5 cm high

Anodized set
(Rear wheel rims & hubs if LIGHT 24",
sideguard supports, forks and fittings) ㉑

- Black
 - Silver
 - Light blue
 - Red
 - Orange
 - Cyclamen light purple
- (only 24" Light wheels)

Options

- 21 "M" velcro strip (hook) on seat
- 22 "F" velcro strip (loop) on seat
- 23 Padded rigid seat h 5
- 24 Chest harness
- 25 Chest belt
- 26 Abductor belt
- 27 Catheter holder pocket
- 28 Rx Lx Crutch holder
- 29 Pump for pneumatics
- 30 Suitcase holder
- 31 45° pelvic belt
- 32 45° padded pelvic belt
- 33 Black holdall
- 34 Umbrella holder
- 35 Front protective covers

Notes
