

COSTA

Hand-bike leggera e versatile, è dotata di un cambio Shimano a 27 rapporti adatto per qualsiasi percorso. La seduta alta, con gambe a 90°, la rende unica per facilità e rapidità di trasferimento ed ideale per praticare ciclismo e divertirsi nel tempo libero.

CARATTERISTICHE TECNICHE

- Telaio in alluminio.
- Ruote da 26" montate con copertoncini.
- Sedile e schienale regolabili in altezza, inclinazione e profondità.
- Pedana regolabile in più posizioni.
- Campanatura ruote posteriori 12° o 6°.
- Comandi dei freni e del cambio direttamente sulle manopole di spinta.
- Doppio freno anteriore.
- Freno di stazionamento.
- Tela schienale tensionabile.
- Fascia reggi gambe e cinturino fermapièdi.
- Paraurti posteriore estraibile (optional).
- Ciclocomputer senza filo.
- Peso: 16 kg (in assetto standard).

COSTA 8CR

Cambio Revo sul mozzo a 8 rapporti interni e freno contropedale.

COSTA 8F

Cambio sul mozzo a 8 rapporti interni con freno a leva.

COSTA 27

Gruppo cambio Shimano a 27 rapporti con guarnitura a 3 corone da 22-28-32 denti, abbinata ad una cassetta a 9 pignoni: serie 11-12-14-16-18-21-24-28-32 oppure serie 12-13-14-15-16-17-19-21-23.

Light and versatile hand-bike fitted with Shimano 27 gears suitable for all terrain. The high sitting position with legs at 90° makes it extremely easy and quick to transfer from/to the wheelchair and ideal for cycling or having fun.

FEATURES:

- Aluminium frame.
- 26" wheels assembled with pneumatic tyres.
- Height, depth, tilt adjustable seat and back-rest.
- Foot-rest adjustable in several positions.
- Rear wheel cambers: 12° or 6°.
- Brake and shift controls on handles.
- Twin front brake.
- Parking brake.
- Width adjustable back-rest upholstery.
- Leg keeper and foot-holding belts.
- Rear detachable bumper bar (optional).
- Wireless computer.
- Weight: 16 kg (standard set-up).



COSTA 8CR

Revo 8 gear hub shift and coaster brake.

COSTA 8F

8 gear hub shift with lever brake.

COSTA 27

Shimano 27 gear shift group with 22-32-44 chainwheel teeth crankset, supplied with a 9 cassette sprockets: 11-12-14-16-18-21-24-28-32 or 12-13-14-15-16-17-19-21-23 sets.



Bauletto portaoggetti.
Top box.



Le foto possono non corrispondere ad ausili in assetto base. / The photos could not show standard version products.